

topspin

**DAISY**

Version 2.0.0



Copyright (C) 2007 by Bruker BioSpin GmbH

revision state: 03/07/2007 (Version 2.0.0)

The information in this manual may be altered without notice.

BRUKER accepts no liability for any mistakes contained in the manual, leading to coincidental damage, whether during installation or operation of the instrument. Unauthorized reproduction of manual contents, in any form, or by any means without written permission of the publishers, or translation into another language, either in full or in part, is forbidden.

Product names used are trademarks or registered trademarks of their respective holders.

Bruker software support is available via phone, fax, e-mail or Internet. Please contact your local office, or directly:

Address: Bruker BioSpin GmbH  
Service & Support Department  
Silberstreifen  
76287 Rheinstetten  
Germany  
Phone: +49 (721) 5161 296  
Fax: +49 (721) 5161 91296  
E-mail: [nmr-software-support@bruker-biospin.de](mailto:nmr-software-support@bruker-biospin.de)  
FTP: <ftp.bruker.de> / <ftp.bruker.com>  
WWW: [www.bruker-biospin.com](http://www.bruker-biospin.com)

---

# Contents

---

<b>Chapter 1</b>	<b>Introduction</b>	<b>5</b>
1.1	About this manual	5
1.2	Conventions	5
1.3	Getting started	6
1.3.1	Basic program features	6
1.3.2	How to start DAISY	6
1.3.3	How to use the DAISY toolbar	7
<b>Chapter 2</b>	<b>DAISY Simulations</b>	<b>11</b>
2.1	Parameters for a DAISY Simulations	11
2.1.1	Open Spin System	12
2.1.2	Importing data	13
2.1.3	The 'Edit Spin System' window	14
2.1.3.1	The Frequency table (of the 'Edit Spin system' window)	14
2.1.3.2	The Scalars table (of the 'Edit Spin system' window)	17
2.1.3.3	The Lineshapes table (of the 'Edit Spin system' window)	19
2.1.4	The 'Advanced Options' window	21
2.1.4.1	The Main table of the 'Advanced Options' window	21
2.1.4.2	The Output Options table of 'Advanced Options' window	22
2.2	Starting a DAISY Simulation	24
2.2.1	Run Simulation/Run Iteration	25
2.2.2	Run Simulation Sequence	25
2.2.3	Generate Lineshape	26
2.2.4	Generate Subspectra	26
2.3	Special Simulation Features	26
<b>Chapter 3</b>	<b>Daisy Iterations</b>	<b>29</b>
3.1	Parameters for a DAISY Iterations	29
3.1.1	Additional parameters in the 'Edit Spin System' window	29
3.1.2	Additional parameters in the 'Advanced Options' window	33
3.1.2.1	Smoothing Parameters	35
3.1.2.2	Additional Options and other Iteration Parameters	36
3.1.2.3	Additional output options in the 'Output Options' window	37

---

<b>Chapter 4</b>	<b>Examples</b>	<b>39</b>
4.1	Spectrum Simulation	39
4.1.1	The <sup>1</sup> H-NMR Spectrum of Aspirine (aromatic protons)	40
4.1.2	The <sup>1</sup> H-NMR Spectrum of Allyl-glycidyl-ether	48
4.2	Spectrum iteration	57
4.2.1	Iteration of the <sup>1</sup> H-NMR Spectrum of Aspirine (aromatic protons)	57
4.2.2	Iteration of the epoxy group of Allyl-glycidyl-ether	64
4.2.3	Iteration of the <sup>1</sup> H-spectrum of trans-Crotonic-acid-ethylester	70
<b>Chapter 5</b>	<b>Appendix</b>	<b>75</b>
5.1	Spin Symmetry, Symmetry groups	75
5.2	Theoretical Background	80
5.2.1	The DAISY Simulator	80
5.2.1.1	The basic Single Spin simulation algorithm	80
5.2.1.2	Using Symmetry	84
5.2.1.3	The Composite Particle Approach	85
5.2.1.4	The nuclei specific linewidth	86
5.2.1.5	The X-Approximation: Division into Subsystems	87
5.2.2	The DAISY Iterator	88
5.2.2.1	The Standard Algorithm	91
5.2.2.2	The Advanced Algorithm	92
5.3	References	93

---

# Chapter 1

## Introduction

---

### 1.1 About this manual

---

This manual is a reference to TOPSPIN structure analysis with 1D spectrum simulation (DAISY). DAISY is a program system designed for simulation and iteration of high resolution NMR spectra. The complete DAISY interface and functionality is described and some examples of simulation will illustrate the power of DAISY.

### 1.2 Conventions

---

#### **Fond conventions**

***daisy*** - commands to be entered on the command line are in Courier bold italic

***Apply*** - commands to be clicked are in Times New Roman bold italic

`fid` - file names are in Courier

contents - any contents of a text file is in Courier small

*name* - any name which is not a filename is in Times New Roman italic

---

## 1.3 Getting started

---

### 1.3.1 Basic program features

The **Spectrum Simulation** is based on the solution of the time - independent Schroedinger equation to calculate the energy levels. To built up the necessary Hamilton matrix the factorisation by total spin values, **chemical equivalence** as well as **magnetic equivalence** and **fragment factorizing** are taken into account to reduce calculation time and memory requirements. Using the selection rules for single quantum transition, frequencies and intensities are determined. Up to 10 nuclei or groups with spin quantum value  $I \geq 1/2$  can be treated with one fragment. The **Lorentzian lineshape** is calculated by application of global or individual line-widths to the transitions of as many fragments as you have defined. The calculated spectrum will be stored in the pdata directory under procno 999 by default.

The **Spectrum Iteration** uses the method of **total-lineshape fitting**. Adjustable parameters are **resonance frequencies**, **scalar coupling constants** and **line-widths**. The program uses the digital experimental spectrum as experimental data.

### 1.3.2 How to start DAISY


As DAISY is now part of **TopSpin**, it will be started as a special analysis option out of the **TopSpin** main window.

When experimental data (processed 1D NMR data) are loaded to the 1D data window, there are two possibilities to start DAISY:

- enter **daisy** into the command line
  - click Analysis -> Structure Analysis -> 1D Spectrum Simulation [daisy] from the menu
- bar

The 1D data window now switches to the DAISY data window.

Figure 1.1: The DAISY data window



### 1.3.3 How to use the DAISY toolbar

The DAISY toolbar contains 22 buttons for input, data preparation, simulation/iteration, display options, scaling, mathematic operation and data handling.



open Spin-System

This button opens a window for selecting files containing spin system information. The default directory is the TopSpin home directory with the default file type setting MGS file (.mgs). (Chapter 2.1.1)



### **import multiplets from a dataset**

With this button the user opens a window where NMR data sets with multiplet can be chosen and imported to DAISY. (Chapter 2.1.2)



### **edit Spin-System**

This button opens a window consisting of three tables: Parameters like shifts/frequencies, couplings, linewidth etc. will be entered here. (Chapter 2.1.3)



### **advanced options**

Pressing this button opens a window again consisting of three tables. Here simulation sequence options, spectrum options, iteration options and output options can be selected. All possible calculations (simulation, iteration, lineshape or subspectra generation) can be started from this window. (Chapter 2.1.4 , 2.2 and 3)



### **run simulation**

Pressing this button will execute a spectrum simulation resulting in a 1D NMR spectrum with the displayed simulated name. (Chapter 2.2.)



### **run iteration**

This command will execute a spectrum iteration. The iteration statistics are available from the protocol list and the resulting 1D NMR spectrum is written. (Chapter 2.2 and 3)



### **show DAISY logfile**

Pressing this button will display the Simulation/Iteration Protocol List



(extension 'lst'/'rst') with the defined windows editor. (Chapter 4.1.1 and 4.2.1., 4.2.2., 4.2.3)



### **deselect all datasets**

This command will deselect all datasets displayed in the DAISY window.



### **reset individual scaling and shifts**

This command resets individual scalings and shifts for the selected dataset.



### **remove all selected datasets from display**

Pressing this button will remove all selected datasets from display.



### **toggle the display layout**

This button switches the display layout from separated spectra to overlaid spectra and vice versa.



### **switch on/off the display of dataset name and scalings**

Pressing this button the dataset name and scaling will be switched on/off.



### **increase intensity**

Pressing this button the intensity of the selected spectra will be doubled.



### **decrease intensity**

Pressing this button the intensity of the selected spectra will be halved.



### **expand intensity axis**

Decreases or increases intensity of selected spectrum while pressing this button and moving the mouse.

**shift up/down**

Shifts up and down the selected spectrum while moving the mouse and pressing this button.

**shift left/right**

Shifts left and right the selected spectrum while moving the mouse and pressing this button.

**show a difference between the 1<sup>st</sup> and the sum of other spectra**

Pressing this button calculates the difference between the first and the sum of other spectra and shows the result as a new spectrum.

**show a sum of all spectra**

Pressing this button calculates the sum of all spectra and shows the result as a new spectrum.

**export Spin-System**

The spin system of the actually loaded DAISY calculation will be stored in the directory <home>NMRSIM/ham (default path) or another directory of your choice. Two files are created (\*.msg, \*.itr)

**terminate DAISY mode, save changes****terminate DAISY mode, do not save changes**

---

# Chapter 2

## DAISY Simulations

---

DAISY calculations can be either simulations or iterations. A simulation means quantum mechanical synthesis of a high resolution NMR spectrum based on given NMR parameters, as there are Resonance frequencies, Scalar Coupling Constants, and linewidths. A spectrum iteration needs all these NMR parameters as well as a starting guess for the iterative refinement of these NMR parameters. In addition experimental data is required for adjustment of the calculated spectrum. As a consequence, both calculation types need the same basic data input, while for iterative calculations simulation data can be extended.

### 2.1 Parameters for a DAISY Simulations

---

This chapter will teach the user to prepare a new DAISY simulation input document by explaining all available features and options, to execute the calculation and display the result.

First of all you have to load the spectrum you want to simulate. It is recommended to choose a special region of your spectrum to simulate. This can be done by using the *int* command. Otherwise the whole spectrum will be used as basic data for a following iteration. This can cause bad results. You should also do peak picking and multiplet definitions using the commands *pp* and *mana*. Then switch to the DAISY window by typing *daisy*.

Now the TopSpin data window looks like it is shown in Figure 1.1. The chosen region of the spectrum is shown on the screen and the toolbar of the data window indicates the DAISY mode.

### 2.1.1 Open Spin System




If simulation data already exists or has been stored with the  button, it can be loaded by clicking the  button on the toolbar. The window shown in figure 2.1 will appear. It shows the TopSpin home directory with the default file type setting MGS File (.mgs).

Figure 2.1: The open-spin-system window



In general DAISY allows to use the following spin system definition formats:

- DAISY file format (\*.mgs file)
- NMR-SIM format (\*.ham file)
- ACD prediction format (\*.acd file)
- Perch prediction format (\*.prp files)

### 2.1.2 Importing data

If you have done a multiplet definition you can import the multiplets from the data set. Multiplet analysis is described in the manual Structure Analysis Tools . DAISY reads out the information from the file `daisymultiplet.txt` in the **procno** directory. Multiplet definition scans the spectrum for identical line distances. It is possible to do this automatically or manually. The more complex the system is, the more manual analysis has to be done. For really complex second order spin systems a normal spin analysis of the spin system is recommended.




Click on the  button in the toolbar and the window shown in figure 2.2 will appear. Here you can type in the file informations or use the browse button. Click the **OK** button. Now you can decide if you want to subdivide the multiplets and if the found coupling constants shall be used.

Figure 2.2: The DAISY import multiplet window



### 2.1.3 The 'Edit Spin System' window

In this window all relevant parameters for the simulation of the chosen spectrum have to be set. Click the  button in the toolbar and a new window will appear. This 'Edit Spin System' window is divided into two parts. The left part is a listing of all fragments for your simulation. It depends on the imported multiplet data, how many fragments will appear. Here fragments can be added, copied or deleted.

The right part consists of three tables, described in the paragraph below:


#### 2.1.3.1 The Frequency table (of the 'Edit Spin system' window)

The 1<sup>st</sup> table (shown in figure 2.3) contains all information about the frequency range of your simulation. It is divided into three blocks.

##### A) Fragment Options:

In the first line a suitable title can be selected. After pressing the **Apply** button this title will be displayed in the listing at the left side of the window.

Figure 2.3: Frequency view of 'Edit Spin System' window



For the creation of the first fragment in a new document you don't need to care about the **Statistical weight** (second line) yet. This value will start to affect the data when more than one fragment is defined. For more information see chapter 2.3.

In the third line a symmetry group can be selected. The default setting for the Symmetry is the point group C1, actually meaning no symmetry at all. As opposed to magnetic equivalence, molecular symmetry leads to chemical equivalence. It is not always necessary to use the maximum symmetry defined by the molecular structure. In general the lower symmetries are much easier to describe and even a simple twofold symmetry greatly improves calculation speed.

Symmetry is useful and saves computer time and memory. Inform yourself about symmetry and gain experimental experience by doing (more in chapter 5.1).

If the selected point group remains C1, the Symmetry Description (available by clicking the symmetry button) can be omitted, because C1 consists of only one symmetry operator, the identity operation. Consequently, further user definitions are not required. If you change the point group selection, all following parameters for this fragment will be affected.

B) Reduce number of lines by:

The default limits of transitions are in normal cases big enough to cover all important lines. The **Minimum intensity** defines the quantum mechanical intensity limit to include a transition into the transition list or not.

C) ISO Value/PSE/Spin Value/Frequencies/Disable:

Below block A and B a frequency block is displayed. It contains the values of the actual fragment. If you have imported data, it is displayed here. Of course, the above declared symmetry directly influences the input.


The list also contains ISO Values, Spin Values, Spins in group and Group Indices. In figure 4 the default values are displayed.

In order to define the nucleus type press the **PSE** button (see figure 4) for each shift/frequency. The 'Choose PSE' window will appear (figure 2.4).

The values from the frequency block are the initial values. Change the element by

pressing the corresponding button in the Periodic table. In the **Nuclear Isotop** combo box the available spinactive nuclei can be selected. According to this selection the element data in the upper right part of the window are updated. Exiting the window by pressing the *ok* button will take the following data to the corresponding nucleus of your spin system: **ISO Value, Element Symbol (on PSE button) and Spin Value**

Figure 2.4: 'Choose PSE' window



The parameters for **Element Symbol** and **Spin Value** are only changeable via the **PSE** button. The **ISO Value** can be changed by user request to apply X-Approximation for homonuclear spin systems, too. X-Approximation will be executed if the element symbol or the ISO value (or both) are different.

The selection of the Isotope and the number of magnetically equivalent spins in a group defines if a single spin or composite particle simulation shall be executed.



The **Resonance Frequency [Hz]/ Chemical Shift [ppm]** values are imported at the start of DAISY or have to be given by the user manually. This is for example necessary when a manual spin analysis has been done. The sum of the spins in group defines the total size of your spin system in the actual fragment.

If we take, for example, an ethoxy-group  $A_2B_3$ , we have two spins in group for A and three spins in group for B.

You can also activate some or all nuclei for iteration (discussed in chapter 4.1.1) and set the upper and lower limits for the frequencies here.

The **Disable** option in the frequency block offers the user to neglect the transitions, which belong to the disabled spin/group. **It does not simulate a decoupling effect!** The coupling information to the disabled nucleus remains in the spin system, only the transitions of the disabled nucleus are suppressed.

This is necessary if you define a fragment in which you have to include one spin/group because of its coupling information towards other spins of this fragment, but for this certain spin/group not all coupling partners are defined. Therefore you want to neglect the transitions of this this special spin/group. This may be possible by restricting the frequency limits. But if the resonance frequency of this special spin/group is located in the same spectral region as the other transitions, it is impossible to suppress its transitions by changing the frequency limits. You have to use the disable options for this spin/group.


Another use for this option is to look at a multiplet pattern of a complex spin system if its transitions are overlapped by other signals.

Due to the fact, that this transition suppression is based on the detection of the transition origin nucleus via eigenvectors, this option may not be exact for strongly coupled systems. Combination lines may not be possible to assign.

In this block you can activate the iterate option for resonance frequencies in the chosen fragment by pressing the iterate *All* button (iteration of all frequencies) or tagging the frequencies you want to iterate. Now it is necessary to set the lower and upper frequency limits. By pressing the *Apply* button all changed values will be stored.

### 2.1.3.2 The Scalars table (of the 'Edit Spin system' window)

The 2<sup>nd</sup> table (shown in figure 2.5) contains all information about the scalar couplings.

Figure 2.5: Scalar Couplings table of 'Edit Spin System' window

This table shows the possible indirect (=scalar) coupling constants. If you have imported the couplings by starting **DAISY** they will appear now. Otherwise you have to fill in your desired values. Correspondingly the selected symmetry affects the coupling table. For fragments with symmetry the minimum number of parameters will be presented. In an unsymmetrical case the maximum number of coupling constants is given.

The coupling constant values are always in Hertz. The table is not influenced by the number of magnetically equivalent spins within a group, because the inner-group scalar coupling constant does not affect the spectrum at all if we consider isotropic solutions only. **DAISY** is not able to calculate NEMA spectra at moment and hence does not read in dipolar coupling constants  $D_{ij}$ . Therefore the quantum mechanical calculation ignores this value and no input is required.

After setting all values press the **Apply** button.

Similar to the frequencies table, you can check iteration of the coupling constants. Again it is necessary then to set lower and upper limits. When choosing a sequence simulation you have to set the number of steps and the step width additionally.

To take over the new values press the *Apply* button.

### 2.1.3.3 The Lineshapes table (of the 'Edit Spin system' window)

The 3<sup>rd</sup> table (shown in figure 2.6) displays all lineshape parameters. By default a **global linewidth** of 0.3 Hz is set.


The **global linewidth** means that this parameter is applied to every transition of the fragment resulting in the Lorentzian lineshape spectrum. The option to use the **Fast Lorentzian Lineshape** is faster but in case of poor digital resolution the resulting spectrum might be inexact. This way the lineshape determination does not calculate every digital point of the spectrum, but picks up values from a model line with the global linewidth. Therefore the transition frequency is moved to the nearest digital point of the spectrum. For this reason the **Fast Lorentzian Lineshape** calculation is not suitable for spectrum iteration.

When checking the option **Nuclei specific linewidth**, the **Fast Lorentzian Lineshape** is no longer available. The display then disables the **global linewidth** and offers to set an **individual linewidth** parameter for every spin/group. The algorithm detects the contribution of the nuclei to every transition based on the eigenvectors.

There are many possible cases where nuclei specific linewidth parameter are needed:

- an hydroxyproton exchanging with water residue
- neighbor nuclei with quadrupole momentum
- unresolved long range couplings

Pressing the *Apply* button will store these parameters.

Figure 2.6: Lineshape Parameters of 'Edit Spin System' window

The input of linewidth has to be given in Hertz. Only independent linewidth parameters are offered for input.

You can check the iteration option like in tables 1 and 2. Now the default iteration limits are available and can be individually modified here.

Press the apply buttons to store the new values.

By pressing the **OK** button you leave the 'Edit Spin System' window.


After setting your individual parameters a simulation can be started by pressing the



button in the toolbar.

## 2.1.4 The 'Advanced Options' window


For more detailed simulations (and, especially, better iterations) you should set some parameters in the 'Advanced Options' window. This window also consists of

three tables. To open this window press the  button in the toolbar.

### 2.1.4.1 The Main table of the 'Advanced Options' window

The first table contains the main options for sequence simulation, iteration options and spectrum options. It is shown in figure 2.7.

Figure 2.7: Main table of the 'Advanced Options' window



The values displayed in figure 2.7 show the default parameters. For a normal simulation you can use these values.


If you enter a frequency offset this value is applied to all resonance frequencies in all fragments when the **Apply** button is pressed.

For more information regarding the Iteration Options see chapter 3.1.2.

#### 2.1.4.2 The Output Options table of 'Advanced Options' window

The second table consists of selectable output options. Figure 2.8 shows the table.

Figure 2.8: Output Options of the 'Advanced Options' window



In this table you can activate simulation and iteration output options as well as result options. In this chapter only the simulation output options will be described. More information about the iteration output option can be found in chapter 3.1.2.3.

#### Output of symmetry data:

**DAISY** can use symmetry factorization (chemical equivalence factorizing) to block out the Hamiltonian matrix more efficiently in order to speed up calculation. Select this item to inspect the character table of the point group symmetry. In case of additional magnetic equivalence (e.g. methyl groups) in your spin system the program lists the symmetry (subs symmetry) of the subspectra.

#### Output of subspectra

This option writes the irreducible representation character of the subspectra into the output protocol in case of symmetry. If Composite Particles are present in the spin system the actual term with weighting of the subspectra is reported.

#### Size of submatrices

In the absence of magnetic equivalence the program writes a table containing the submatrices for every irreducible representation arranged according to descending total spinvalues. The last column shows the total spin factorization in case of no selected symmetry. If magnetic equivalence is present the program detects the size of submatrices at runtime, so the information is given in the output protocol, when the appropriate matrix is composed (not at the beginning of the calculation).

#### Linear combinations

This option gives a total representation of the symmetrized basic spin functions used for setup of the Hamiltonian matrix. These linear combinations are NOT the eigen vectors, they are just the symmetrical functional basis.

#### Eigenvalues

All calculated eigenvalues are saved into the protocol file in the order of their determination.

### Transitions and Energy levels

For every submatrix the transition frequency and intensity and the pair of corresponding energy levels are listed.

### Frequencies before degeneracy

Towards the end of the calculation all frequencies obeying the quantum mechanical restrictions as defined are listed in decreasing order. In addition a file with the extension 'fre' is created containing the same list. The existence of this file enables the additional **Run** command **Generate Lineshape**. Executing this command will not lead to a quantum mechanical recalculation of the spectrum, but the 'fre'-file containing the determined transition list will be read in. The actual global line width of the corresponding fragment will be used to create the simulated NMR-spectrum. Be careful, only the lines of a fragment for which this option is checked, is listed in the file.

### Frequencies after degeneracy

Similar to the latter option, but the transitions belonging to the same frequency are combined and corresponding intensities are summed up. No 'fre' file is written.

In the '**Result block**' of this table you can set the Destination PROCNO. It is recommended to choose 'Increment'. The simulation gets the next possible PROCNO of the experimental spectrum. If you check 'Fixed' you can choose a PROCNO. This will be used for every simulation if you don't change the value.

---

## 2.2 Starting a DAISY Simulation


---

After setting the output options, switch to the 3<sup>rd</sup> table, the 'Run table' of the 'Advanced Options' window, which is shown in figure 2.9. Press the **Run Simulation** button in this table or leave the 'Advanced options' window by using the **OK**

button and press the  button in the toolbar.





Figure 2.9: 'Run table' of the 'Advanced Options' window



The 'Run table' offers five possibilities to calculate spectra with DAISY.

### 2.2.1 Run Simulation/Run Iteration

The buttons *Run Simulation* and *Run Iteration* correspond to the  - and  - buttons in the toolbar of the DAISY window.

### 2.2.2 Run Simulation Sequence

The *Run Simulation Sequence* button allows to calculate spectra with systemati-

cally changed NMR parameters (resonance frequencies and/or coupling constants). The number of steps and the step width for frequencies must have been activated by pressing the *Seq. Sim.* button and checking the frequencies in the Frequencies table (see Chapter 2.1.3.1.). This calculation is called Run Simulation, because quantum mechanical calculations will take place for every individual calculation. The first simulation will be the reference calculation - without variation of any parameter.

### 2.2.3 Generate Lineshape

If you have selected the Output Option **Frequencies before degeneracy** (see Chapter 3.1.4.2.) a list of transitions and intensities of this fragment is written to a 'fre'- file. After a simulation/iteration with this option checked, pressing the **Generate Lineshape** button reads this list and generates a lineshape based on the defined spectrum parameters and the global line-width of the lineshape parameters of this fragment. Even if the option 'nuclei specific linewidth' has been selected, this command cannot deal with specific lineshape parameters, because no eigenvectors (on which the individual linewidths are detected) are available. This command does not take care of if you have selected the output option for all the fragments. It just takes the transitions of these fragments, for which the option has been selected. For small spin systems this calculation is not faster than a quantum mechanical calculation, because the lines have to be read from the harddisk. For large fragments, however, (e.g. 10 spins) this option might be useful.

### 2.2.4 Generate Subspectra

In case the spin system contains either Composite Particles (groups of magnetically equivalent nuclei and/or Spin Values  $>1/2$ ) and/or Symmetry (for definition of symmetry refer to section (2.1.3.1)).

## 2.3 Special Simulation Features

---

There are different situations conceivable where the option *statistical weight* is indispensable for spectrum simulation. It is only effective if more than one fragment is defined:

Calculation of an unequal mixture of compounds in the solution, e.g. isomeric mixtures, impurity signals.

When a compound contains not only pure, spin-active elements, the spectrum consists of the main and a satellite spectrum. The latter one includes the not 100%-abundant spin interactions.

If parts of the molecule give equal patterns they may be defined just once and weighted according to their presence.



---

# Chapter 3

## Daisy Iterations

---

DAISY iteration input documents are based on simulation input files with some additional information required. The basic features and definition of indispensable parameters will be described. References to the DAISY simulation will be mentioned if the features are already described there.

### **3.1 Parameters for a DAISY Iterations**

---

The following chapter will teach you setup and start a DAISY iteration. If you have simulated the spectrum you can use the input parameters of this simulation. Some additional parameters for an iteration have to be added.

If you don't have simulation parameters, start DAISY and follow the parameter setup as described in chapter 2.

#### **3.1.1 Additional parameters in the 'Edit Spin System' window**

For iteration setup some additional inputs in all the tables of the 'Edit Spin System' window are required.

In the 'Frequencies' table, the spin(s) you want to iterate have to be checked. Now

it is possible to fill in the lower and upper limit for the chemical shift or resonance frequency (depending on what is checked ppm or Hz). An example for such a table is shown in figure 3.1.

Figure 3.1: 'Frequencies' Table of 'Edit Spin System' window with iteration parameters

Fragment Options

Fragment Title: 2x AMX System mit ch2-koppl.  
 Statistical weight: 1.0000  
 Symmetry: C1

Reduce number of lines by  
 Lower limit for transitions: -1999.48 [ppm]  
 Upper limit for transitions: 1999.48 [ppm]  
 Minimum intensity: 0.0010

Add/	Def	Disable	ISO	PSE	Spin	Iterate	Lower limit	Chemical shift	Upper limit	Spins in	Group	Annotation
			Value	Value	Value	None	[ppm]	[ppm]	[ppm]	group	Index	
X	<input type="checkbox"/>	1	H	H	1/2	<input checked="" type="checkbox"/>	3.7463	3.7544	3.7663	1		
X	<input type="checkbox"/>	2	H	H	1/2	<input checked="" type="checkbox"/>	3.4104	3.4307	3.7749	1		
X	<input type="checkbox"/>	3	H	H	1/2	<input checked="" type="checkbox"/>	3.1713	3.1914	3.2113	1		
X	<input type="checkbox"/>	4	H	H	1/2	<input checked="" type="checkbox"/>	2.8099	2.8298	2.8499	1		
X	<input type="checkbox"/>	5	H	H	1/2	<input checked="" type="checkbox"/>	2.6262	2.6463	2.6661	1		

Seq. Sim. >>>

Switch for iterate All or None

Checkbox for activating an iteration

If iteration is activated you can fill in here the frequency/shift limits

PPM C:HZ

OK Apply Cancel

- **Group Index**

Some NMR spectroscopists may not want to define chemical equivalence via Point Group Symmetrie or due to other reasons parameters shall be kept identical for optimization although they are not tied together by symmetry. Therefore

a feature has been added to DAISY allowing to tie parameters of the same type together, with the purpose to keep them identical. Valid types are Resonance Frequencies and Coupling Constants.

**IT IS NOT POSSIBLE TO COMBINE SYMMETRY AND GROUPING!**

To visualize the difference, if parameters are dependent either because of chemical equivalence or because of manual grouping, in the latter case the edit fields are shown, but disabled. In case of symmetrical dependent parameters these fields are invisible. For grouped parameters only the first appearing group member is enabled for editing, all others are disabled. Their connection can be seen from the identical number in the Group Index entry. The number zero means no grouping.

The group index is a nonzero integer number and is valid for the actual parameter type over all fragments. This means, parameters can be grouped together over all fragments of a DAISY document. For example all resonance frequencies having the group index '1' in the whole DAISY document will get the same properties, except the disable option in the 'Frequencies table'. Properties are:

- parameter values
- iteration options
- iteration limits

and for frequencies the same number of:


- spins per group
- spin value
- isotope value
- element symbol etc.

Group indices can be set in the 'Frequencies' table and/or the 'Scalars Couplings' table.

In the 'Scalar Couplings' table the iteration of the coupling(s) must be activated, too. This can be done by checking the box beneath 'Iterate'. Now the fields with lower and upper limits are available and the desired values can be filled in. As distinct from the 'frequencies' table with its check box for each frequency, here you have to choose the coupling first by clicking on it. The value of the chosen coupling is now the 'Start Value' of the iteration parameters. Checking the iteration

box activates this coupling for iteration. An example is shown in figure 3.2 .

Figure 3.2: 'Scalars Couplings' Table of 'Edit Spin System' window with iteration parameter




In the 'Lineshapes' table you have the possibility to do Lorentzian Lineshape calculation with a global linewidth or a calculation with a Nuclei Specific Linewidth (see Chapter 2.1.3.3). An iteration with Lorentzian Lineshape calculations iterates all lines of the spectrum with the same global linewidth. If you check the **Nuclei Specific Linewidth** box, an iteration for all or only some spins with their specific



linewidth is available. An example is shown in figure 3.3. **Do not use Fast Lorentzian Lineshape for spectrum iteration!**

Figure 3.3: 'Lineshapes' Table of 'Edit Spin System' window with iteration parameter



For each table the new value will be taken over by pressing the **Apply** button. To leave this window press then the **OK** button.


### 3.1.2 Additional parameters in the 'Advanced Options' window

The three tables of the 'Advanced Options' window contain some iteration specific parameters, too. They are explained in the following section.

A DAISY iteration is more time efficient when iterating only regions containing NMR signals of the spins of interest. To define the regions for the signal groups, each signal group in your experimental spectrum should be integrated and the integration saved.

For a DAISY iteration it is possible to import these integration regions in the 'Main' table of the 'Advanced Options' window. Press the **Iter.Regions** button (see figure 2.7 in chapter 2.1.4.1). A new popup window will appear as shown in figure 3.4.

Figure 3.4: Window with import data for iteration regions



Select the regions for the iteration and press **Import Regions** button.

In the 'Iteration options' block of the 'Main' table you can select your desired iteration mode. Two modes are available, both of them are total-lineshape fitting methods.

- **Standard**

The Standard iteration mode represents the original spiral algorithm written by Binsch et al, which was optimized by Boenigk<sup>9d</sup>, based on the combination of

steepest descent and Gauss-Newton corrections for the parameters. (Chapter 5.2.2.1)

- **Advanced**

The Advanced modus uses successive steepest descent and Gauss-Newton correction, optionally with a double-sum-target function, introduced by Höffken<sup>11</sup>. (Chapter 5.2.2.2)

By default the Standard iteration mode is set.

### 3.1.2.1 Smoothing Parameters

The parameters **Broadening** and **No. of Cycles** in this options block are connected to the Iteration mode. The combination of these parameters define the smoothing of the spectra. A brief explanation of the meaning of these parameters without going deep into theory:

Total-lineshape fitting is designed to avoid letting the iteration fall into local minima and to reach the global minimum. The starting parameters for iteration can be quite poor. The poorer the parameters the stronger the starting smoothing has to be in order to enable the system to find correct solution. E.g. if your simulated spectrum with your initial parameter guess for iteration is already quite good (or even almost correct), use only low smoothing (or even none).

- **Broadening**

Possible values are **none**, **low**, **medium** and **high**. **None** means, the original experimental data will directly be used for iteration without smoothing.

**Low**, **medium** and **high** define the initial strength of spectrum smoothing.

The stronger the starting smoothing is selected, and the more points the spectral regions contain, the longer the iteration will take in the first iteration cyclus.

- **Number of Cycles**

You can select **low**, **medium** and **high**.

This value defines the factor for reducing the smoothing during iterations. Here the same effect is valid: The higher the value is selected, the longer the iteration will take, but the worser the starting NMR parameters can be.

### 3.1.2.2 Additional Options and other Iteration Parameters

- **Iterate baseline**

Checking this option will include the correction of the experimental spectrum with a baseline increment and a baseline tilt in the last phase of iteration. But this baseline iteration parameter does not substitute a proper baseline correction with the TopSpin processing command **bas** before using **DAISY**.

In case of the Standard Iteration Mode an additional parameter may be selected:

- In case the relative differences of resonance frequencies are well known, it is useless to iterate all resonance frequencies independently. For this purpose the iteration of only one resonance frequency (the first in the first fragment) shall be selected for iteration and the option **Iterate Frequency offset** shall be checked. When this option is used ONLY the very first resonance frequency in the first fragment must be selected for iteration, NO other frequency is allowed to be selected. During optimisation the frequency differences from the very first shift to all other shifts is kept constant! This feature is e.g. very useful to analyse spectra of compounds in anisotropic solutions, when isotropic NMR parameters are well known.

In the case of the **Advanced Iteration Mode** an additional parameter may be selected:

- **Starting Frequencies poor**

If this option is checked, the Advanced iteration will start with the double-sum target function in the first phase.

There are two other parameters which are important for iteration:

- **Convergence Criterion**

The convergence criterion (default 0.05, means the difference between consecutive error values) defines the limit, when to change to the next iteration cycle to decide no further improvements can be made with the actual smoothing.

Normally this value doesn't need to be changed.

- **Default number of iterations**

This value can be entered to stop the iteration, when no minimum can be reached. For normal cases this number of iterations won't be reached.

### 3.1.2.3 Additional output options in the 'Output Options' window

In the 'Output Options' window there are five additional output options besides the simulation output options.

- **Correction vectors**

For every iteration the parameter vector together with the gradient vector and (Standard)/ or (Advanced) the Gauss-Newton vector is combined and written into the output file.

- **Search for error reduction**

For Standard iteration three different spirals with four testing points on each spiral are used and for Advanced iteration the distance is halved to achieve error reduction. This option informs about the spiral geometry and the actual error for each test point.

- **Store intermediate results**

Especially when a time-consuming iteration is started, this option is very useful. When it is activated the resulting parameters of every single iteration are stored in a file with a consecutive number as part of the file extension. So instead of file name <name>.mgs, the first iteration result will be <name>.mgs1, the second <name>.mgs2 etc.. Be careful, because 100 iterations will generate 100 files in the procno directory of the experimental NMR spectrum. If you check this option you have to delete these files manually, if you don't need them anymore.

- **Covariance matrix**

The parameter covariance matrix is written in the output file.

- **Correlation matrix**

The correlation matrix is written in the output file.

Now all parameters are set and you can start an iteration by pressing the







button in the toolbar or out of the 'Run' table in the 'Advanced Options' window. In the procno-directory of NMR-spectrum (<TopSpin>\pdata\procno) 4 files will be created by >DAISY:

- daisy.itr
- daisy.log

- daisy.lst
- daisy.mgs

If you want to keep these files, you have to rename them. DAISY will overwrite them when starting a new simulation/iteration of the chosen spectrum.

After you have done your DAISY calculations it is possible to store the calculated spin system by pressing the  button in the toolbar of the DAISY window. By default two files (\*.msg, \*.itr) will be stored in TopSin home directory. Any other directory can also be chosen, too. These files contain the data of the tables in the 'Edit Spin System' window. The stored data can be loaded by pressing the  button when needed.

To terminate the DAISY session press the  button or the  button. The first one terminates DAISY saving changes, the second one terminates without saving changes.

---

# Chapter 4

## Examples

---

This chapter illustrates the calculation of experimental spectra with DAISY. A suitable setup will be shown for each example.

Section 4.1 will illustrate the DAISY simulation, whereas section 4.2 will explain some features of parameter optimization by DAISY iteration.

### 4.1 Spectrum Simulation

---


For all simulations and/or iterations you have to treat your experimental spectrum in the same way before entering DAISY.

- After normal processing (ft, pk, cal) you need a peak picking for the whole spectrum or for the signal groups you want to simulate (see TopSpin Users Guide). Without this the multiplet definition will not work.
- It is recommended to integrate the different signal groups of your spectrum. You get information about the number of protons in a signal group. For a possible iteration after your simulation you will need the integration to define your iteration regions (see TopSpin Users Guide).

### 4.1.1 The $^1\text{H}$ -NMR Spectrum of Aspirine (aromatic protones)

The first example of a DAISY simulation is the aromatic proton system of aspirine. The structure is shown in figure 4.1.

Figure 4.1: Structure of aspirine




The  $^1\text{H}$ -NMR spectrum shows four separated signal groups representing the four aromatic protons of the compound. In order to prepare the spectrum for a DAISY simulation a peak picking is requested in the first step. Afterwards it is possible to execute a multiplet analysis with Multiplet Definition (*mana*). As the signal groups are clearly separated and show good resolution, the automated multiplet analysis is fine. There is no need for manual corrections. The automatic multiplet definition starts at the low-field side of the spectrum. The numbering of the multiplets was adapted to the protons in the structure using the information from other NMR experiments (see figure 4.2.).


Using the multiplet report available in the Multiplet Definition window you can assign connections between the found shifts and couplings. (see figure 4.3 ). These determined connections will be transferred to the 'scalars' table in the 'Edit Spin System' window when starting DAISY.



Figure 4.2.:  $^1\text{H}$ -Spectrum of aromatic protons of aspirine with peak picking and multiplet definition




For identical couplings found in several multiplets, the average of these values will be filled in. Looking to the list in figure 4.2, you will realize that there are no couplings  $J(1,4)$  and  $J(2,3)$  visible. The coupling from H<sub>1</sub> to H<sub>4</sub> is a long range coupling over 5 atoms. It can be set to 0. The coupling  $J(2,3)$  is an ortho-coupling like  $J(2,4)$  and  $J(1,3)$ . It must be of the same size like these couplings. This missing value will be set manually in the 'scalars' table of the 'Edit Spin System' window. Before leaving the multiplet definition window save all data by clicking the 'terminate/save changes' button.

Figure 4.3. Assign Connection/Report in from Multiplet Definition window

After starting DAISY you will be asked if you want to import information from the multiplet definition. Click Ok and the following data will appear in the "Edit Spin System window":

Figure 4.4: 'Frequency' table for the aspirine aromatic protons



As multiplet definition had found 4 multiplets, a 4-spin system is built. We clear the iterate check boxes for a first simulation. They are always checked when starting with complete data (frequencies and couplings) from multiplet definition.

Figure 4.5: Scalar Couplings table for aspirine aromatic protons





Figure 4.6: Lineshape table for aspirine aromatic protons





After importing all data, the coupling  $J(2,3)$  has to be filled in manually (see above). For a first simulation a global line with of 0.3 Hz is ok. Start the simulation by pressing the  button. Figure 4.7 shows the simulated and measured spectrum.

Figure 4.7: Simulated and measured spectrum of aspirine aromatic protons



Looking at the result we see both resonance frequencies and couplings fit quite well. But the linewidth is not correct. In the simulation the lines of the different multiplets are much better resolved than in the real spectrum. Besides this we need line specific linewidth. The input data of this simulation are quite well for a final iteration of this system (see chapter 4.2.1).

The DAISY report of this simulation contains input data and data chosen in the output options. In this example it was 'Frequencies before and after degeneration'.

The following paragraphs show the DAISY log-file which appears pressing the



button in the DAISY window. You can find this file named `daisy.lst` in the `/procno-directory` of your spectrum.

```
SPIN SYSTEM CONSISTS OF 4 NUCLEI
// SINGLE SPIN SIMULATION
=====
```

Standard Output:  
detected by the programm

```
TITLE: Fragment 1: 4-spins in 4 group(s)
-----
STARTPARAMETERS DEFINING CALCULATED SPECTRUM
=====
```

Standard output:  
the fragment title  
Standard Output

```
ISO-VALUE  RESONANCE FREQUENCY
-----
1          4083.5864
1          3835.0218
1          3700.9119
1          3592.5588
```

Standard Output:  
complete resonance frequencies

```
COUPLINGS  SCALAR  DIPOLAR
-----
1 2        1.5504   0.0000
1 3        7.7520   0.0000
1 4        0.0000   0.0000
2 3        7.9000   0.0000
2 4        7.9021   0.0000
3 4        1.5500   0.0000
```

Standard Output:  
complete scalar couplings

```
GLOBAL LINEWIDTH : 0.3000
-----
```

Expected Number of Calculated Transitions: 56

Standard output:  
the total intensity of all possible  
and of all calculated lines is given here

```
LINES CALCULATED:          56  32.000
```

LINES SUPPRESSED BY

-----

- MINIMUM INTENSITY ( 0.001): 24 0.000  
 - LOWER FREQ. LIMIT (-1000000.000): 0 0.000  
 - UPPER FREQ. LIMIT ( 1000000.000): 0 0.000

LINES BEFORE DEGENERACY: 32 32.000

**Frequencies before degeneracy**  
 Transition frequencies and intensities  
 ordered for decreasing frequencies

1 seconds

0.00

FREQUENCY INTENSITY LINENUMBER

-----

4088.2795	0.9737	29
4088.2784	0.9738	26
4086.7228	0.9858	20
4086.7213	0.9857	13
4080.5344	1.0136	23
4080.5339	1.0137	15
4078.9777	1.0268	11
4078.9767	1.0267	4
3843.8750	0.9156	30
3842.3183	0.9079	17
3835.9820	0.9739	27
3835.9750	1.0321	24
3834.4249	0.9671	9
3834.4183	1.0143	7
3828.0821	1.1035	16
3826.5248	1.0857	3
3709.3637	1.0631	31
3707.8049	1.0999	28
3701.6187	1.0172	18
3701.4638	0.9495	21
3700.0604	1.0554	10
3699.9050	0.9685	14
3693.7186	0.9128	5
3692.1604	0.9337	2
3597.2153	1.0476	32
3597.2142	1.0486	19
3595.6565	1.0173	25
3595.6559	1.0165	8

3589.3223	0.9803	22
3589.3207	0.9809	6
3587.7635	0.9547	12
3587.7625	0.9539	1

LINES AFTER DEGENERACY:

=====

**Frequencies after degeneracy**  
intensities of degenerated lines are added  
up

FREQUENCY	INTENSITY	NUMBER
-----		

4088.2795	0.9737	1 +
4088.2784	0.9738	2 +
4086.7228	0.9858	3 +
4086.7213	0.9857	4 +
4080.5344	2.0274	5 +
4078.9777	1.0268	6 +
4078.9767	1.0267	7 +
3843.8750	0.9156	8 +
3842.3183	0.9079	9 +
3835.9820	0.9739	10 +
3835.9750	1.0321	11 +
3834.4249	0.9671	12 +
3834.4183	1.0143	13 +
3828.0821	1.1035	14 +
3826.5248	1.0857	15 +
3709.3637	1.0631	16 +
3707.8049	1.0999	17 +
3701.6187	1.0172	18 +
3701.4638	0.9495	19 +
3700.0604	1.0554	20 +
3699.9050	0.9685	21 +
3693.7186	0.9128	22 +
3692.1604	0.9337	23 +
3597.2153	1.0476	24 +
3597.2142	1.0486	25 +
3595.6565	2.0339	26 +
3589.3223	0.9803	27 +
3589.3207	0.9809	28 +
3587.7635	0.9547	29 +
3587.7625	0.9539	30 +

LINES AFTER DEGENERACY: 30


1 seconds

0.01

### 4.1.2 The $^1\text{H}$ -NMR Spectrum of Allyl-glycidyl-ether

The second example for DAISY simulations is one of the two clearly separated 5-Spin-systems of Allyl-glycidyl-ether. The structure of this compound is shown in figure 4.8.

Figure 4.8: Structure of Allyl-glycidyl-ether



Between 2.5 ppm and 3.8 ppm there are five separated multiplets which represent the three protons of the epoxy-ring and the CH<sub>2</sub>-group.

For a DAISY simulation chemical shift/resonance frequency informations and the coupling constants of the protons are needed. The analysis option **Multiplet Definition** can help to find these parameters. For the protons 1 - 4 (see figure 14) this can be done by automatic multiplet definition. Each of the four protons shows a dublett of dubletts, so there is no problem to find couplings and chemical shifts.

For proton 5 the situation is more difficult. It is in a center position of this group and therefore shows couplings to all four protons. The linewidth of the signal group for proton 5 indicates that some lines are probably not resolved. The multiplet definition is difficult and may lead to wrong couplings and chemical shifts.


The multiplet definition of the other groups contains all required couplings. So we only read out the chemical shift for proton 5 (center of the signal group) and leave the 'Multiplet Definition' window while storing the results in the file `multi-plet.txt`.

Now we type *daisy* in the command line or choose Analysis/Structure Analysis/1D Spectrum Simulation from the main menu bar of TopSpin. Figure 4.9 shows



the actual display.


Figure 4.9: Import of data for DAISY calculations (1<sup>st</sup> step)



Click **OK** in the popup window and DAISY will import the multiplet information you created with 'Multiplet Definition'.

After clicking **OK** (or **Cancel**) button a new popup window will appear. Now you have to decide where you want to subdivide your system into independent fragments. (see figure 4.10)


In our example we cancel this to calculate our 5-spin-system as one fragment.

Figure 4.10: 2<sup>nd</sup> step of data import for DAISY calculation

Finally you can to import unassigned couplings from the `multiplet.txt` file. (see figure 4.11)

We don't do it for this example because we have to add a frequency and the corresponding couplings manually as explained above.

Figure 4.11: 3rd step of data import for DAISY calculation



Now all data gets imported to DAISY. We can check and complete the frequencies and couplings in the different tables of the 'Edit Spin System' window.


Press the  button in the toolbar of the Data window and the frequencies table will appear (see figure 4.12). The first four frequencies/shifts have been imported from the file `multiplet.txt`, the fifth frequency has been set manually.

Figure 4.12: Frequencies table for the epoxy-CH<sub>2</sub>-group of Allyl-glycidyl-ether

the four imported frequencies


manually added frequency

press Apply button to store data

Del	Disable	Value	PSE	Spin	Iterate	Lower limit	Chemical shift	Upper limit	Spins in	Group
X		Γ 1	H	1/2	Γ	1.7393	3.7560	1.7103	1	
X		Γ 2	H	1/2	Γ	1.4293	3.4295	1.4425	1	
X		Γ 3	H	1/2	Γ	1.0133	2.0350	1.0007	1	
X		Γ 4	H	1/2	Γ	1.0263	2.6465	1.0000	1	
+		Γ 5	H	1/2	Γ	1.1723	3.1921	1.1121	1	


Afterwards we have to check and complete the Scalar coupling table (see figure 4.12). As we did not import unassigned couplings, only coupling  $J_{12}$  is set in the table. The other couplings are entered in manually.

Figure 4.12: Scalar Couplings table for the epoxy-CH<sub>2</sub>-group of Ally-glycidyl-ether



The last table, the Lineshape table is important for our example. As explained above, proton 5 (see figure 4.8) shows different couplings to the four other protons of the group. The resulting signals are closely together and don't seem to be resolved. Therefore the signal group representing proton 5 shows linebroadening. For the simulation of the whole spin system we use a global linewidth of 0.5 Hz.


Before starting the simulation some output options and the Destination PROCNO are set. Figure 4.13 shows the chosen values.

Figure 4.13: Output Options for the epoxy-CH<sub>2</sub>-group of Allyl-glycidyl-ether


The simulation of the spectrum using these parameters is shown in figure 4.14. Although it is 'only' a first simulation, experimental and simulated spectrum fit together very well. A closer look at the signal group for proton H<sub>5</sub> (figure 4.15) shows a well chosen linewidth of 0.5 Hz.

These simulation parameters are good starting values for an iteration of this 5-Spin-System. They are stored together with your chosen output options in the simulation protocol in the file `daisy.lst`.

Figure 4.14: Experimental and simulated spectrum of the epoxy-CH<sub>2</sub>-group of Allyl-glycidyl-ether



Resonance frequencies and couplings seem to fit quite well. The global linewidth of 0.5 Hz seems to be better justified for Allyl-glycidyl-ether than for the previous case (aspirine) is better than the one in the first example (aspirine). nevertheless an advanced calculation using spin specific linewidth is requested. (see chapter 4.2.2).

Figure 4.15: Experimental and simulated spectrum of proton H<sub>5</sub>

The frequencies, scalar couplings and linewidth used for the simulation again can be found in the log-file, the `daisy.lst`, respectively in the /procno-directory of your spectrum.



## 4.2 Spectrum iteration

The DAISY iteration input documents are based on simulation input files with some additional information.

### 4.2.1 Iteration of the $^1\text{H}$ -NMR Spectrum of Aspirine (aromatic protons)

The 4-spin system of the aromatic protons of aspirine shows four clearly separated multiplets in the  $^1\text{H}$ -NMR-spectrum. Multiplet analysis and DAISY simulation with determined frequencies and couplings are described in chapter 4.1.1. The simulation results showed that the frequencies and couplings are quite well. So we can use these values as start input for an iteration. Looking at the simulation we realized the need of a nucleus specific linewidth for an exact calculation of this spectrum.

Figure 4.16: Frequency table prepared for iteration

The screenshot shows the 'Frequency View' window with the following data:


Add/	Del	Disable	ISO	Spin	Iterate	Lower limit	Chemical shift	Upper limit	Spins in	Group	Annotation
			Value	PSE	Value	[ppm]	[ppm]	[ppm]	group	Index	
X	<input type="checkbox"/>	1	H	1/2	<input checked="" type="checkbox"/>	8.1451	8.1651	8.1850	1		
X	<input type="checkbox"/>	2	H	1/2	<input checked="" type="checkbox"/>	7.6481	7.6681	7.6880	1		
X	<input type="checkbox"/>	3	H	1/2	<input checked="" type="checkbox"/>	7.3799	7.3999	7.4199	1		
X	<input type="checkbox"/>	4	H	1/2	<input checked="" type="checkbox"/>	7.1633	7.1633	7.2032	1		

A red arrow points to the 'Iterate' column, which is checked for all four protons, indicating that their frequencies will be iterated.

For the following iteration all parameters available in the 'Edit Spin System' window will be iterated. Therefore click on the iteration check boxes for each frequency, scalar coupling and linewidth. The figures 4.16, 4.17 and 4.18 show the three tables of the 'Edit Spin System' window.

The upper and lower limits for iteration are automatically set by DAISY to +10 Hz and -10 Hz to the start value (see figure 4.16). These values can be changed manually. Using fine start values, this is not necessary.

Figure 4.17: Scalar coupling table prepared for iteration




Similar to the frequencies, the upper and lower limits for iteration of the couplings are set automatically by DAISY. For proper chosen start values these limits are fine. In our example the long range coupling 1,4 is not resolved in the experimen-

tal spectrum. So we will leave this value at 0 Hz and do not iterate it.


In the lineshape table we choose Nuclei Specific Linewidth and iterate all nuclei. Again DAISY automatically sets the upper and lower iteration limits.

Figure 4.18: Lineshape table prepared for iteration




In the 'Advanced Options' window the additional parameters have to be set as described in chapter 3.1.2. To get good results we import the iteration regions of the experimental spectrum and check 'Iterate BaseLine' (Main table of this window). Figure 4.19 shows the upcoming window for definition of the iteration regions.

Figure 4.19 Importing iteration regions




'Start' and 'End' in this table are the start and end points of the integration region of the experimental spectrum. The most relevant information for the spectrum shape is available in the signal regions. So it is only necessary to iterate the spectrum in these regions. The baseline in between will be calculated. Now the iteration can be started. When it is done you will be asked if you want to transfer the calculated values in the input tables (see figure 4.20).


Figure 4.20: DAISY window after end of iteration



Clicking **OK** transfers the calculated values and frequencies etc. to the tables of the 'Edit Spin System'. Now a new iteration could be started with these new start values. Otherwise the input values will be kept. The calculated values can be seen in the log-file.

Figure 4.21: Iterated spectrum overlaid with experimental spectrum



The following paragraph shows the start parameters and the result of the iteration. It was taken from the file `daisy.rst` in the directory `<user>/pdata/1`. This is also the source for the data available with the  button in the DAISY window.

## Parameters to be iterated

```

=====
Fragment  Type   starting value  lower limit  upper limit
-----
1  F( 1)      4083.59000    4073.58640   4093.58640
1  F( 2)      3835.02000    3825.02180   3845.02180
1  F( 3)      3700.91000    3690.91190   3710.91190
1  F( 4)      3592.55000    3582.55880   3602.55880
1  J( 1, 2)    1.55040      -0.44960     3.55040
1  J( 1, 3)    7.75200       5.75200     9.75200
1  J( 2, 3)    7.92100       5.90000     9.90000
1  J( 2, 4)    7.92100       5.90000     9.90000
1  J( 3, 4)    1.00000      -0.45000     3.55000
1  H( 1)       0.30000       0.00100     2.00000
1  H( 2)       0.30000       0.00100     2.00000
1  H( 3)       0.30000       0.00100     2.00000
1  H( 4)       0.30000       0.00100     2.00000
0  Baseline Offset  0.00000     -2.00000     2.00000
0  Baseline Ascent  0.00000     -2.00000     2.00000

```

Iteration 21 rms = 1.78270396

Index Parameter Type Parameter Value neg. Grad. GN-Corr.

```

-----
0  F( 1)      4083.65011    2.891e-004   2.952e-004
1  F( 2)      3834.98820   -2.802e-004  -2.614e-004
2  F( 3)      3700.91634   3.834e-004   3.713e-004
3  F( 4)      3592.53425   3.880e-004   3.802e-004
4  J( 1, 2)    1.69580       7.011e-004   1.503e-003
5  J( 1, 3)    7.87206      -6.621e-004  -3.382e-004
6  J( 2, 3)    7.39957       2.575e-003   1.005e-003
7  J( 2, 4)    8.09810      -5.682e-004   1.634e-004
8  J( 3, 4)    1.08665      -9.393e-004   2.475e-004
9  H( 1)       1.54867      -4.263e-003  -4.733e-003
10 H( 2)       1.34961      -5.573e-003  -5.830e-003
11 H( 3)       1.29116      -3.844e-003  -3.774e-003
12 H( 4)       1.71547      -7.271e-004  -1.068e-003

```

## 4.2.2 Iteration of the epoxy group of Allyl-glycidyl-ether

As mentioned in chapter 4.1.2, the simulation input parameters lead to a simulated spectrum which fits well to the experimental one (see figure 4.14 and 4.15). A final iteration of the parameter will lead to a well simulated NMR-spectrum. In this example, the definition of lower and upper limits is of primary importance to obtain a successful results. The values that are automatically by DAISY define a wide range of parameter variation. It is too big for this example with its 5-spin system and the very different linewidth for each nucleus. But if you confine the range for the parameters, especially for the frequencies and the linewidth, the iteration will be fast and fine. Otherwise the iteration leads to a deterioration of the simulated spectrum. The system probably falls into a local minimum. (see chapter 5.2.2.). The figures 4.22 - 4.24 show the input tables for a successful iteration.

Figure 4.22: Frequencies table prepared for iteration of epoxy-CH<sub>2</sub>-group of Allyl-glycidyl-ether

Add	Del	Disable	ISO	Value	PSE	Spin	Iterate	Lower limit [ppm]	Chemical shift [ppm]	Upper limit [ppm]	Spins in group	Annotation
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	Γ 1	1	H	1/2	F	3.7500	3.7500	3.7600	1	
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	Γ 2	1	H	1/2	F	3.4255	3.4295	3.4395	1	
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	Γ 3	1	H	1/2	F	2.8267	2.8308	2.8347	1	
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	Γ 4	1	H	1/2	F	2.6429	2.6468	2.6509	1	
<input type="checkbox"/>	<input checked="" type="checkbox"/>	<input type="checkbox"/>	Γ 5	1	H	1/2	F	3.1881	3.1921	3.1961	1	

The upper and lower limits of the frequencies, scalar couplings and specific linewidth are set manually and preferably close to the start parameters. This is possible because the simulation shows that the input values are nearly the true values. Otherwise the simulated spectrum would not fit that well to the experimental one.



Figure 4.23: Scalar coupling table prepared for iteration of epoxy-CH<sub>2</sub>-group of Allyl-glycidyl-ether





Figure 4.24: Lineshape table prepared for iteration of epoxy-CH<sub>2</sub>-group of Allyl-glycidyl-ether




With these parameters and the iteration regions shown in figure 4.25 we start a DAISY iteration.

Figure 4.25: Iteration regions window (available under 'Advanced Options')



Start	End
6.0056	5.8680
5.3846	5.1565
4.1650	3.9644
3.8272	3.6871
3.4945	3.3518
3.2418	3.1376
2.8794	2.7968
2.6897	2.5962

In the iteration regions window the iteration regions of the whole experimental spectrum are available because the whole spectrum was integrated at the beginning of the analysis (see chapter 4.1.2). We choose the five regions between 3.8 and 2.5 ppm and press the button 'Import Regions'. Then the iteration can be started by

pressing the  button in the toolbar of the DAISY window or out of the 'Run' table of the 'Advanced Options' window. Figure 4.26 shows the result. To illustrate the quality of the iteration the figure also shows the experimental spectrum.

As mentioned above the signal groups in the experimental spectrum have different linewidths. Especially the signal group representing Proton H<sub>5</sub> shows a significant line broadening compared to the other signals. Iterating with specific linewidths for each nuclei leads to a satisfying result which can be seen more detailed in figure 4.27.

Figure 4.26: Experimental and iterated spectrum of epoxy-CH<sub>2</sub>-group of Allyl-glycidyl-ether






Figure 4.27: Experimental and iterated signals of H<sub>5</sub>-proton of epoxy-group of Allyl-glycidyl-ether



The detailed results of iteration can be seen clicking the  button in the menu bar or opening the file `daisy.rst` in `<user>/pdata/1`. Input values and final iteration value are shown now:

Iteration

=====

Maximal number of iterations = 100  
 Convergence criterion = 0.05000  
 Correlation factor = 0.01000

Multiplier = 6.00000  
 Total number of cycles = 4  
 Correlation factor lineshape = 0.10000  
 Iteration of lineshape starting in the 3rd Cycle

Parameters to be iterated

=====

Fragment	Type	starting value	lower limit	upper limit
1	F (1)	1877.98820	1875.49000	1880.49000
1	F (2)	1715.19580	1713.20000	1717.20000
1	F (3)	1415.76800	1413.74000	1417.74000
1	F (4)	1323.74410	1321.77000	1325.77000
1	F (5)	1596.46500	1594.46000	1598.47000
1	J (1, 2)	-11.45000	-12.45000	-10.45000
1	J (1, 5)	3.05000	2.05000	4.05000
1	J (2, 5)	5.85000	4.85000	6.85000
1	J (3, 4)	-4.90130	-5.95000	-3.95000
1	J (3, 5)	4.25110	3.25110	5.25110
1	J (4, 5)	2.70000	1.70000	3.70000
1	H (1)	0.50000	0.00100	1.20000
1	H (2)	0.50000	0.00100	1.20000
1	H (3)	0.50000	0.00100	1.20000
1	H (4)	0.50000	0.00100	1.20000
1	H (5)	0.50000	0.00100	1.20000
0	Baseline Offset	0.00000	-2.00000	2.00000
0	Baseline Ascent	0.00000	-2.00000	2.00000

Statistical Information :

=====

Final sum of squares = 297.999358  
 Number of spectral points = 1622  
 Degrees of Freedom = 1606  
 Standard Deviation of Measurements = 0.430760  
 R-Factor ( %) = 3.918803

Fragm. No.	Parameter type	Initial Parameter Vector	Best Parameter Vector	Standard Deviation
1	F (1)	1877.98820	1878.26984	0.00072

1	F( 2)	1715.19580	1715.30507	0.00072
1	F( 3)	1415.76800	1415.66754	0.00099
1	F( 4)	1323.74410	1323.81176	0.00096
1	F( 5)	1596.46500	1596.55872	0.00243
1	J( 1, 2)	-11.45000	-11.45124	0.00101
1	J( 1, 5)	3.05000	3.07087	0.00137
1	J( 2, 5)	5.85000	5.86312	0.00135
1	J( 3, 4)	-4.90130	-5.05153	0.00127
1	J( 3, 5)	4.25110	4.17002	0.00154
1	J( 4, 5)	2.70000	2.73139	0.00176
1	H( 1)	0.50000	0.78791	0.00142
1	H( 2)	0.50000	0.78866	0.00141
1	H( 3)	0.50000	0.93455	0.00198
1	H( 4)	0.50000	0.94549	0.00191
1	H( 5)	0.50000	1.07432	0.00489


2 seconds

0.83

### 4.2.3 Iteration of the $^1\text{H}$ -spectrum of trans-Crotonic-acid-ethylester

The last example that will be analysed with DAISY is the proton spectrum of the trans-Crotonic-acid-ethylester (s.figure 4.28.).

Figure 4.28.: Structure of Crotonacidethylester



The  $^1\text{H}$ -spectrum of this substance shows 5 clearly separated and very well dissolved signal groups (H2 - H6, see annotation in DAISY frequencies table, too). According to the structure of the molecule it is a good example for dividing the spectrum in two independent subspectra. The Ethylester group is totally independent from the system around the double bond. A coupling of the protons via carboxy

group is not possible.

Figure 4.29: Frequencies table for iteration of Fragment 1 (Vinylgroup)






Figure 4.30: Frequencies table for iteration of Fragment 2 (Vinylgroup)




For both fragments a corresponding 'Scalar Couplings' table and a 'Lineshape' table is available. The values are filled in directly from the multiplet definition file `daisy multiplet.txt`. For the following DAISY calculation of the spec-

trum all parameters (frequencies, couplings and linewidth) are checked for iteration. In the 'Lineshape' table **Nuclei Specific Linewidth** is chosen and iterated for each spin. The result is shown in figure 4.31.

Figure 4.31 Simulated and original spectrum of trans-Crotonic-acid-ethylester



Again the detailed results of iteration can be seen clicking the  button in the menu bar or opening the file `daisy.rst` in `<user>/pdata/1`. Input values and final iteration value are shown now:



## Iteration

=====

Maximal number of iterations = 100

Convergence criterion = 0.05000

Correlation factor = 0.00010

Multiplier = 6.00000

Total number of cycles = 7

Correlation factor lineshape = 0.02154

Iteration of lineshape starting in the 4th Cycle

## Parameters to be iterated

=====

Fragment	Type	starting value	lower limit	upper limit
1	F(1)	2755.16510	2745.16510	2765.16510
1	F(2)	2303.48840	2293.48840	2313.48840
1	F(3)	717.24300	707.24300	727.24300
2	F(1)	1637.97220	1627.97220	1647.97220
2	F(2)	478.85560	468.85560	488.85560
1	J(1, 2)	15.49500	13.49500	17.49500
1	J(1, 3)	6.95560	4.95560	8.95560
1	J(2, 3)	1.72390	-0.27610	3.72390
2	J(1, 2)	7.16570	5.16570	9.16570
1	H(1)	0.30000	0.00100	2.00000
1	H(2)	0.30000	0.00100	2.00000
1	H(3)	0.30000	0.00100	2.00000
2	H(1)	0.30000	0.00100	2.00000
2	H(2)	0.30000	0.00100	2.00000

## Statistical Information :

=====

Final sum of squares = 109.752807

Number of spectral points = 1505

Degrees of Freedom = 1491

Standard Deviation of Measurements = 0.271312

R-Factor ( % ) = 1.062955

Fragm. Parameter type Initial Parameter Best Parameter Standard

---

No.	Vector	Vector	Deviation	
1	F (1)	2755.16510	2755.06463	0.00279
1	F (2)	2303.48840	2303.65185	0.00322
1	F (3)	717.24300	717.15556	0.00104
2	F (1)	1637.97220	1637.93539	0.00115
2	F (2)	478.85560	478.80272	0.00083
1	J (1, 2)	15.49500	15.55628	0.00400
1	J (1, 3)	6.95560	7.01139	0.00181
1	J (2, 3)	1.72390	1.87399	0.00159
2	J (1, 2)	7.16570	7.18705	0.00110
1	H (1)	0.30000	1.10007	0.00547
1	H (2)	0.30000	1.14072	0.00634
1	H (3)	0.30000	1.33525	0.00209
2	H (1)	0.30000	1.25424	0.00230
2	H (2)	0.30000	1.41031	0.00167

2 seconds

1.03

---

# Chapter 5

## Appendix

---

This chapter gives some more information about symmetry groups and symmetry description (chapter 5.1) and a short summary about the theoretical background of the DAISY simulator and the DAISY iterator (chapter 5.2)

### 5.1 Spin Symmetry, Symmetry groups

---

The spin symmetry in DAISY is described using the Schoenflies permutation group. The default input is  $C_1$  which means that the molecule does not contain any symmetric element. For the non-symmetry  $C_1$  no input of symmetry description is required. All other symmetries need permutation operators to fix the molecular symmetry. The Schoenflies symbol have been chosen here, because they are often used in the description of molecular symmetry.

To give an overview which symmetry requires how many permutations, the following table is given. The first row gives the total number of permutation operators including the unity operation, so that one permutation less than enumerated has to be defined.

number of permutations	Schoenflies point group symmetries
1	$C_1$
2	$C_2, C_s, C_i$
3	$C_3$
4	$S_4, D_2, C_4, C_{2v}, C_{2h}$
5	$C_5$
6	$S_6, D_3, C_6, C_{3v}, C_{3h}$
8	$D_{2d}, D_{2h}, C_{4h}, C_{4v}, D_4$
10	$C_{5h}, C_{5v}$
12	$T, D_{3d}, D_{3h}, D_6, C_{6v}, C_{6h}$
16	$D_{4d}, D_{4h}$
20	$D_{5d}, D_{5h}$
24	$O, T_h, T_d, D_{6h}, D_{6d}$
48	$O_h$
60	$I$

**Table 5.1** : Schoenflies point groups and corresponding permutations

The Symmetry Description dialog box directly depends on the Symmetry Group selection made. For default  $C_1$  Symmetry Group no userdefined Symmetry Description is necessary.


The Symmetry Description dialog box informs about the number of descriptions, which means the number of permutation operators. This value minus 1 represents the number of lines required for permutation input. This defines the layout of the dialog box.

In the headline the point group name and the number of corresponding symmetry description is displayed. The unity-operation is given in the first line and is at the

same time the fixed nucleus numbering. In the first row the name of the symmetry operation is displayed to the user preparing the input properly. It is not allowed to interchange the sequence of symmetry operations. This will lead to errors during calculation and wrong results.

For symmetry higher than C1 all permutations have to be defined before closing the dialog box. This is necessary because the construction of the following input tables (for frequencies/couplings) depends on the symmetry definition.

The creation of permutation operators is described using naphthalene as example:



The spin system of naphthalene can be either described as  $[AB]_4$  ( $C_{2v}$ ) or as  $[AB]_4$  ( $D_{2h}$ ). We will show the creation of the permutation table with the symmetry group  $C_{2v}$  and explain that this is the better choice, although  $D_{2h}$  strikes the maximum molecular symmetry.

The  $C_{2v}$  symmetry consists of 4 symmetry operators: the identity (E), the  $C_2$ -axis ( $C_2$ ) and two vertical symmetry planes (SigV and SigV')

Figure 5.1: Initial Symmetry Description for naphthalene



The identity results in the subsequent spin numbering. The other three symmetry

operators are followed by the initial numbering of the spin system (see figure 5.1).

Now the  $C_2$  operation will be applied to the molecule according to the numbering of the spin system leading to the permutation operators. The  $C_2$  axis is perpendicular to the molecular plane and interchanges spins number 1 into 5, number 2 into 6, number 3 into 7, number 4 into 8 and vice versa (see figure 5.2).

Figure 5.2: Symmetry Description completed with  $C_2$  permutation

The screenshot shows a dialog box titled "Symmetry Description" with a dropdown menu set to "C2v" and "No. Of Descriptions" set to 4. Below this is a table with 4 rows and 8 columns. The columns are labeled 1 through 8, and the rows are labeled E, C2, SigV, and SigV'. The C2 row shows a permutation: 1 to 5, 2 to 6, 3 to 7, and 4 to 8.

	1	2	3	4	5	6	7	8
E	1	2	3	4	5	6	7	8
C2	5	6	7	8	1	2	3	4
SigV	1	2	3	4	5	6	7	8
SigV'	1	2	3	4	5	6	7	8

The next operation involves the SigV plane. This (and all subsequent symmetry operations) operate on the spin labelling as defined in the first line (identity operation E). Spins 1 and 4 change places, and so do spins 2 and 3, 5 and 8, and finally 6 and 7 (see figure 5.3).

Figure 5.3: Symmetry Description completed with SigV permutation

The screenshot shows the same dialog box as Figure 5.2, but the SigV row in the table now shows a permutation: 1 to 4, 2 to 3, 3 to 2, 4 to 1, 5 to 8, 6 to 7, 7 to 6, and 8 to 5.

	1	2	3	4	5	6	7	8
E	1	2	3	4	5	6	7	8
C2	5	6	7	8	1	2	3	4
SigV	4	3	2	1	8	7	6	5
SigV'	1	2	3	4	5	6	7	8

The fourth symmetry operation  $\text{SigV}'$  is now perpendicular to the  $\text{SigV}$  plane. The effect of this mirror plane transfers spin 1 to 8, 2 to 7, 3 to 6 and 4 to 5 (see figure 5.4).

Figure 5.4: Symmetry Description completed with  $\text{SigV}'$  permutation

The screenshot shows a dialog box titled "Symmetry Description" with a close button (X). It displays "Symmetry Group" as "C2v" and "No. Of Descriptions" as "4". Below this is a table with 4 rows and 8 columns. The columns are labeled 1 through 8, and the rows are labeled E, C2, SigV, and SigV'. The table contains the following permutation values:

	1	2	3	4	5	6	7	8
E								
C2	5	6	7	8	1	2	3	4
SigV	4	3	2	1	8	7	6	5
SigV'	8	7	6	5	4	3	2	1


At the bottom right of the dialog box are "OK" and "Cancel" buttons.

As mentioned before, the spin system of naphthalene satisfies all required symmetry operations of  $D_{2h}$  (identity, three perpendicular  $C_2$  axes, one inversion center and three perpendicular mirror planes - every plane is constructed of two  $C_2$  axes). So there are twice as many symmetry operators in  $D_{2h}$  than in  $C_{2v}$ .

Looking at the additional permutation information about the symmetry of  $D_{2h}$  in the naphthalene case (see figure 5.1.5), we will recognize that the permutations are in pairs identical ( $E = \text{SigXY}$ ,  $C_{2z} = i$ ,  $C_{2y} = \text{SigYZ}$ ,  $C_{2x} = \text{SigXZ}$ ). This does not yield additional information about the molecule.

For DAISY calculation this does not result in a gain of calculation speed. As this is the main reason for symmetry consideration, there is no need to use  $D_{2h}$  symmetry group instead of  $C_{2v}$  symmetry group in this case.

Conclusion: In many cases it is not necessary to select the highest symmetry that is available in the Symmetry Description window. A partially reduced symmetry will do for fast calculation.

Figure 5.5: D<sub>2h</sub> Symmetry Description for naphthalene


E	1	2	3	4	5	6	7	8
C2Z	5	6	7	8	1	2	3	4
C2Y	4	3	2	1	8	7	6	5
C2X	8	7	6	5	4	3	2	1
i	5	6	7	8	1	2	3	4
SgXY	1	2	3	4	5	6	7	8
SgXZ	8	7	6	5	4	3	2	1
SgYZ	4	3	2	1	8	7	6	5

## 5.2 Theoretical Background

References: 1-18

### 5.2.1 The DAISY Simulator

References: 1-9,12,16-18

#### 5.2.1.1 The basic Single Spin simulation algorithm

References: 8,9,17

The simulation of high resolution NMR spectra is based on the time independent Hamilton Operator (here in the first step for Single Spins) to satisfy the Schrodinger equation in the Hilbert space.



$$\hat{H} = - \sum_{i=1}^N \frac{1}{2\pi} \gamma_i h B_0 (1 - \sigma_i) I_{z,i} + \sum_{i < k}^N \sum_{i < k}^N h J_{ik} \hat{I}_i \hat{I}_k \quad (1)$$

**H** nuclear spin Hamilton operator [Joule]

**N** number of nuclei in spin system

$\gamma_i$  gyromagnetic ratio of the i-th nucleus

**h** Planck constant

**B<sub>0</sub>** static magnetic field (defines in -z direction)

$\sigma_i$  shield constant of the i-th nucleus

**I<sub>z,i</sub>** spin operator of the i-th nucleus in z direction

**J<sub>ik</sub>** Scalar (indirect Coupling Constant between i-th and k-th nucleus)

**I<sub>i</sub>** spin operator of the i-th nucleus

After dividing the Hamiltonian operator by **h**, using the definition of the resonance frequency

$$\nu_i = \frac{1}{2\pi} \gamma_i B_0 (1 - \sigma_i) \quad (2)$$

$$\hat{H} = - \sum_{i=1}^N \nu_i \hat{I}_{z,i}$$

$$\begin{aligned}
 & \sum_{i < k}^N \sum^N (J_{ik}) \hat{I}_{z,i} \hat{I}_{z,k} + \\
 & \frac{1}{2} \sum_{i < k}^N \sum^N (J_{ik}) \hat{I}_i \hat{I}_{-k} + \hat{I}_{-i} \hat{I}_k
 \end{aligned} \tag{3}$$

$H$  nuclear spin Hamilton operator [Hz]

$I_{+,i}$  raising operator of spin  $i$

$I_{-,i}$  lowering operator of spin  $i$

called the Zeeman term (field dependent), the coupling term (including the z-spin operator) and the mixing term (including raising and lowering term).

For first order calculation the X-approximation, which neglects the mixing term of the Hamiltonian, is implemented.

The set of basic functions is composed of the product spin function  $|\Phi_i\rangle$ , defined as the product of the spin functions  $|s_k\rangle$ :

$$|\Phi\rangle = \prod_{k=1}^N |s_k\rangle \tag{4}$$

$N$  = number of nuclei in spin system

For single spin calculations involving spins with  $I=1/2$ , the spin functions  $|sk\rangle$  are symbolised as  $|\alpha\rangle$  or  $|\beta\rangle$ .

The eigenfunctions  $|\Psi_i\rangle$  to solve the Schroedinger equation

$$\hat{H}|\Psi_i\rangle = E_i|\Psi_i\rangle \quad (5)$$

are derived as linear combinations of basic product spin functions.

$$|\Psi_i\rangle = \sum_{j=1}^{2^N} c_{ij}|\Psi_j\rangle \quad (6)$$

The dimension of the Hamilton matrix is given by the number of different basic function  $2^N$  identical with the number of eigenvalues  $E$ .

The eigenvalues  $E$  are formed by diagonalisation of the Spin Hamiltonian Matrix. Diagonalisation of matrices is a very time consuming process. Consequently matrices are broken down to submatrices of smaller dimension thus efficiently reducing the time demand. The standard way of factorisation is to use the commutator

$$[\hat{H}, \hat{F}_z] = 0 \quad (7)$$

with  $F_z$ , the total spin operator, defined by

$$\hat{F}_z |\Phi_k\rangle = m_{T,k} |\Phi_k\rangle \quad (8)$$

to obtain the total spin factorization.

### 5.2.1.2 Using Symmetry

References: 5,8,9,16,17

Further factorisation is achieved by introducing symmetry considerations. A commutator is used:

$$[\hat{H}, \hat{P}_n] = 0 \quad (9)$$

where  $P_n$  stands for the projection operator resulting from the n-th irreducible representation of the spin symmetry group. The basic functions  $|\Phi_i\rangle$  are linear combinations of the product spin functions to create the eigenfunctions:

$$|\tilde{\Phi}_i\rangle = \sum_{j=1}^{2^N} c_{ij} |\Phi_j\rangle \quad (10)$$

The coefficients  $c_{ij}$  are based on the information of the character table of the point group for spin symmetry.

Using the new basic function  $|\psi\rangle$

$$|\Psi_{k'}\rangle = \sum_{j=1}^{2^N} c_{kj} |\tilde{\Phi}_j\rangle \quad (11)$$

the eigenfunctions are derived as required by the Schroedinger equation.

Each Hamiltonian sub-submatrix is diagonalised individually. Then the allowed transition frequencies and intensities are determined according to known selection rules. The demands for memory and processing time largely depend on the size of these sub-submatrices. So symmetry factorisation is important for saving computer time.

Of course the main factor governing the total computer time will be determined by the number of nuclei involved in the spin system. The advanced DAISY user will quickly learn to select the relevant spin-islands to efficiently calculate certain regions of the total spectrum.

### 5.2.1.3 The Composite Particle Approach

References: <sup>8,16</sup>

If magnetical equivalence occurs and/or spins with spin values higher than 1/2 are involved, the Composite Particle approach is automatically used for speedy calculations. Here every Composite Particle group is treated like one particle in every term of the system.

The properties of the Composite Particle will be described by the following concept:

The particle is decomposed into a weighted superposition of individual states, which will be dealt with as a Single Spin of  $I \geq 1/2$  and handled as individual quantum mechanical systems. The resulting subspectra are finally combined to give a total spectrum of the spin system. Note that the quantum mechanics to be calculated for composite particles is nearly the same as for single spins. The only difference is, that the basic functions  $|sk\rangle$  can be composed of any possible spin value, not just only  $|\alpha\rangle$  and  $|\beta\rangle$ .

Symmetry can be applied for the main term - the molecule itself. For all other terms the program detects contained subgroups of symmetry automatically from the implemented subgroup-tree. If symmetry is present some of the terms give identical spectra, so that the spectrum will be calculated once and then summed up including the proper statistical weights as derived from binomial coefficients.

The Hamiltonian for the Composite Particle approach for every subspectrum is of the same type like the Single Spin Hamiltonian of equation (3). The only difference is, that for spin system calculation in anisotropic solution the Hamiltonian can be extended.

#### 5.2.1.4 The nuclei specific linewidth

Reference:<sup>16</sup>

Arising from various effects (neighboring of quadrupolar nuclei, dynamic effects etc.) the presence of different linewidths in one spectrum can often be seen. To solve this problem in the most general way for all possibilities - even for second order systems and overlapping spectral regions - the assignment of specific linewidth to every nucleus is the best theoretical approach.

Based on these specific line width parameter (index i)  $HWB_i$  for every resonance (index j) in the spectrum a line specific width  $hwb_j$  will be determined which can be traced back to the eigenvector coefficients  $c_i$  for the actual transition (j).

$$hwb_j = \sum_{i=1}^N c_i HWB_i \quad (12)$$

For normalization the eigenvectors are re-scaled in the following way:

$$\sum_{i=1}^N |c_i|^2 = 1 \quad (13)$$

### 5.2.1.5 The X-Approximation: Division into Subsystems

Reference: <sup>18</sup>

For increasing size of spin systems the dimension of the Hamiltonian matrix is growing rapidly. Various methods to break them down into smaller parts have been developed over the years:

- the factorisation by the Total-spin values
- the factorisation by permutational operators
- the factorisation by different molecular terms (Composite Particle approach)

The 'simple' X-Approximation based on different ISO-Values solely makes off-diagonal elements of the Hamiltonian matrix to vanish. The dimension of the matrix itself is not touched at all.

The basic equation for the X-Approximation is the perturbation parameter  $\lambda_{ik}$ .

$$\lambda_{ik} = \left| \frac{J_{ik}}{\nu_i - \nu_k} \right| \quad (14)$$

As an approximate value for neglect, the mixing term from the Hamiltonian  $\lambda_{ik}$  should be much smaller than 0.01. For different isotopes the numerator (difference between the Resonance Frequencies) is in the order of  $10^{-6}$  Hz (MHz), so that this condition is always fulfilled for different types of nuclei, even if a Coupling Constant of  $10^3$  Hz is present. For nuclei of the same isotope the limit for X-

Approximation is more difficult to fix.

However a very simple case occurs when the Coupling Constant is zero. Then independent from the denominator the numerator is zero. The absolute value of a spectral parameter is zero by definition and that's the smallest value obtainable:

$$\lim \lambda_{ik} = \lim \left| \frac{0}{v_i - v_k} \right| = 0 \quad (15)$$

But the neglectation of the off-diagonal elements is not the only consequence, also the diagonal elements are not affected by any Coupling Constant contribution. The sole influence of the other nucleus is an additive contribution to all energy levels, so that the transition frequencies are not affected.

An extreme situation occurs when a whole cluster of nuclei present the property against another cluster of nuclei. Here the 'intelligent' X-Approximation starts:

Real factorization according to the membership of a cluster (subsystem, fragment) leads to independent quantum mechanical systems - with much smaller matrix-dimension - in which the former described factorization takes place. To deal with these systems different input cases can be combined to a complete input document. The calculation will be executed in a sequence of simulation runs for all the fragments.

### 5.2.2 The DAISY Iterator

References: <sup>8-18</sup>

DAISY includes a powerful simulation and iteration algorithm. The iteration may be classified as a total line shape fitting or - more generally - as an intensity iteration. This approach is different from traditional frequency iterations based on LAOCOON type programs.

No trial-simulation and no line assignments are required for iterations with DAISY. The iteration environment uses efficient simulators for Single Spins and Composite Particles described in the preceding chapter.

The iteration algorithm is a modification of previous results obtained or used by



Binsch et al.<sup>15)</sup> in DAVINS, e.g. the iteration routine is based on the spiral algorithm published by Jones<sup>14)</sup>. Previous DAISY versions<sup>8,9)</sup> were designed in FORTRAN IV and FORTRAN77 to run on mainframe computers and work stations. The FORTRAN versions were able to handle up to 8 spins. Considerable progress was made with WIN-DAISY, the new modified C-version using totally dynamic memory allocation. The 1994 version covers up to 10 spins/groups (Composite Particles) per fragment. The current DAISY version, now implemented in Top-Spin, is based on this 1994 version of WIN-DAISY, which was developed as a cooperative project between the research team of G. Hägele in Düsseldorf and Bruker Spectrospin, Karlsruhe<sup>8-12,16-18)</sup>.

As mentioned above the iteration routine uses the spiral algorithm, which is based on previous studies from D.W. Marquardt<sup>13)</sup> and K. Levenberg<sup>13)</sup>. The error of the initial parameter guess (which has to be minimized) is calculated using a Gaussian least square treatment for all spectra data points.

$$\sum_{i=1}^n (res_i)^2 \quad (16)$$

n        number of data points  
res<sub>i</sub>    residual in the i-th data point

The residual is defined as:

$$res_i = f_i(p) - exp_i \quad (17)$$

f<sub>i</sub>(p) = calculated intensity in the i-th data point  
exp<sub>i</sub> = experimental intensity in th i-th data point  
p     = NMR parameter

The Jacobian matrix  $J$ , consisting of first differentials of the Lorentzian line shape and the NMR parameters, are determined.

$$J = \begin{bmatrix} \frac{\partial}{\partial p_1} f_1(p) & \frac{\partial}{\partial p_n} f_1(p) \\ \frac{\partial}{\partial p_1} f_k(p) & \frac{\partial}{\partial p_n} f_k(p) \end{bmatrix} \quad (18)$$

$n$  = number of data points

$k$  = number of parameters to be iterated

$p$  = NMR parameter

The Hessian matrix, the second derivative matrix, is built by multiplying an exponential weighting vector to the approximate Hessian, briefly defined by the expression  $J^T J$ .

$$Hess_{appr} = J^T J \quad (19)$$

The Gauss-Newton correction is defined as

$$c_{Gauss-Newton} = -Hess_{appr} J^T res \quad (20)$$

$c_{\text{Gauss-Newton}}$	=	Gauss-newton correction
$\text{Hess}_{\text{appr}}$	=	approximated Hessian matrix
$J^T$	=	transposed Jacobian matrix
res	=	residual

The negative gradient can be defined as:

$$c_{\text{neg, grad}} = -J^T \text{res} \quad (21)$$

The corrections applied to the NMR parameter are composed of negative gradients and Gauss-Newton correction, which can be done in different ways.

In many cases, sharp local minima occur in the total line shape fitting, consequently an efficient smoothing adopted to the iteration algorithm is required which can be done again in different ways.

### 5.2.2.1 The Standard Algorithm

References: 8,9,12,15,16

Iteration theory derives, that one of the best error reductions in these complicated cases is situated in the plane between the Gauss-Newton correction and the negative gradient.

Corrections applied to parameters are thus not simple Gauss-Newton corrections used in traditional frequency iterations, but composed of the Gauss-Newton vector and the negatively scaled gradient.

The error reduction starts with applying the total length of the Gauss-Newton correction to the parameters under iteration. If no reduction of the error is achieved, a spiral is constructed and three more testing points in the plane are inspected to reduce the error value. If this is not successful, the same procedure will be retried with other test spirals starting from a lower Gauss-Newton correction (reductions: division by a factor of three).

Altogether twelve error reduction tests are made, if necessary, for every iteration. The effective correction can be defined as:

$$c_{eff} = \lambda c_{negGrad} + (1 - \lambda) c_{Gauss - Newton} \quad (22)$$

with  $0 \leq \lambda \leq 1$

The differentials forming the Jacobian matrix are calculated numerically. For each parameter two spectral simulations with small changes in the actual parameter values in positive and negative direction are performed. This iteration procedure is valid for the linewidth as well. Base line increment and base line tilt are determined directly.

### 5.2.2.2 The Advanced Algorithm

References: <sup>10,11,17</sup>

As it is essential to have properly adjusted main spectral parameters for the use of the Standard Algorithm, a novel double sum target function (17) was developed. This function was created to be used with conjugated gradients ("starting frequencies poor").

$$\sum_{i=1}^n \sum_{j=1}^i (res_j)^2 = min$$

The result of this iteration is not exact, so the standard target function with quasi-Newton minimization is applied afterwards. The difference to the former described Gauss-Newton method is the substitution of the Hessian matrix by proper matrices to avoid the inversion procedure.

This calculation is faster than the standard method but there is no general rule, which approach (Standard or Advanced) would be the optimum to reach the cor-

rect result. Even more, there is no exact solution to predict the iterability of a given spectrum.

### 5.3 References

---

- [1] J.A. Pople, W.G. Schneider and H.J. Bernstein; "High Resolution Nuclear Magnetic Resonance", McGraw Hill, New York, 1959
- [2] P.L. Corio, *Chem. Rev.*, **60** (1960) 363
- [3] P.L. Corio; "Structure of High Resolution NMR Spectra", Academic Press, New York, 1966
- [4] R.K. Harris; "Nuclear Magnetic Resonance Spectroscopy", Pitman, London, 1983
- [5] R.G. Jones; "The Use of Symmetry", in "NMR-Basic Principles and Progress", Vol. 1, p. 97ff, Springer, Berlin-Heidelberg-New York, 1969
- [6] C.W. Haigh; "A Simple Guide of the Use of Iterative Computer Programs in the Analysis for NMR Spectra", in "Annual Reports on NMR Spectroscopy", Vol. 4, Academic Press, London, 1971
- [7] S.M. Castellano, "LAOCN3" in "Computer Programs for Chemistry", Vol. 1, p. 10ff, Benjamin, New York, 1968
- [8] G. Hägele, M. Engelhardt and W. Boenigk, "Simulation und automatisierte Analyse von Kernresonanzspektren", VCH, Weinheim, 1987
- [9] a) V. Lueg and G. Hägele; *J. Mag. Res.*, **26** (1977) 47  
b) V. Lueg; Dissertation, Universitaet Düsseldorf, 1977  
c) W. Boenigk and G. Hägele; *QCPE Bulletin* **4(1)** (1984) 20  
d) W. Boenigk; Dissertation, University Düsseldorf, 1984  
e) G. Hägele, W. Boenigk and M. Engelhardt; *Frezenius Z. f. Anal. Chemie*, **320** (1985)676  
f) G. Hägele, P. Reinemer and M. Grzonka; in "Computer in der Chemie

- Softwareentwicklung in der Chemie" **2** (1988) 241
- g) S. Goudetsidis and G. Hägele; in "Computer in der Chemie Softwareentwicklung in der Chemie " **4** (1990) 233
- [10] B. Neeb; Diploma; University Düsseldorf, 1988
- [11] H.-W. Höffken; Diploma, Heinrich-Heine-University Düsseldorf, 1990
- [12] U. Weber; Diploma, Heinrich-Heine-University Düsseldorf, 1991
- [13] a) D.W. Marquardt, *SIAM J. Appl. Math.*, **11** (1963) 431  
b) K. Levenberg, *Quart. Appl. Math.*, **2** (1944) 164
- [14] A. Jones, *The Computer Journal*, **13** (1970) 301-307
- [15] a) D. Stephenson and G. Binsch, *J. Magn. Reson.*, **37** (1980) 395  
b) D. Stephenson and G. Binsch, *J. Magn. Reson.*, **37** (1980) 395
- [16] U. Weber, Dissertation, Heinrich-Heine-University Düsseldorf, 1994
- [17] H.-W. Höffken, Dissertation, Heinrich-Heine-University Düsseldorf, 1994
- [18] R. Spiske, Dissertation, Heinrich-Heine-University Düsseldorf, 1995

---

# Index

---

## A

- advanced options 8, 21
  - advanced iteration mode 35, 36
  - frequency offset 22
  - iterate frequency offset 36
  - iteration options 21
  - iteration regions 34
  - selectable output options 22
  - sequence simulation 21
  - spectrum options 21
  - standard iteration mode 36
  - standard iterations mode 34

## C

- calculated spectrum
  - adjustment 11
  - iteration 6, 11
  - simulation 6, 11
- composite particles 26

## D

- DAISY data window 7
- DAISY iterator 88
  - advanced algorithm 92
  - standartd algorithm 91
- DAISY log-file of asperine simulation 45
- DAISY Simulator
  - Composite Particle Approach 85
  - Single spin algorithm 80
  - symmetry consideration 84
  - X-Approximation 87

## dataset

- deselect 9
  - remove from display 9
  - switch on/of display of name and scaling 9
- ## destination procno 24, 37

## E

- element symbol 16
- equivalence
  - chemical 6, 15
  - magnetic 6, 15

## F

- Fast Lorentzian Lineshape 19
- fragments
  - listing 14
  - options 14
  - statistical weight 15
  - subdivision of multiplets 13, 49
- frequency table 14, 42, 52

## I

- import multiplets 8, 42, 48
- In 31
- input document 11
  - importing data 13
  - loading 12
- ISO Value 16
- iterate options
  - activate couplings 19, 31
  - activate lineshape parameters 20
  - activate resonance frequencies 17, 29
  - group index 30

iteration of Asperine spectrum  
  experimental spectrum and iteration 61  
iteration input document 29  
  setup of 29  
iteration of Asperine spectrum 57  
  definition of iteration regions 60  
  experimental spectrum and iteration 62  
  frequency table 57  
  lineshape table 59  
  scalar coupling table 58  
  transferring the calculated values to  
    'Edit Spin System' tables 61  
iteration of Crotonacidethylester spectrum 70  
  experimental spectrum and iteration 72  
  frequencies table for the fragments 71  
  iteration results 73  
iteration of epoxy group of Allyl-  
glycidyl-ether spectrum 64  
  definition of iteration regions 66  
  experimental spectrum and iteration 67  
  frequencies table 64  
  iteration results 68  
  lineshape table 65  
  scalar coupling table 65

## **L**

lineshape parameters 20  
  table for aromatic protons of Asirine 43  
linewidth 6  
  global 19  
  nuclei specific 19, 86  
logfile 8  
Lorentian lineshape 6  
Lorentzian lineshape calculation  
  with global linewidth 32  
  with nuclei specific linewidth 32

## **M**

Minimum intensity 15

## **O**

output  
  additional options 37  
  correction vectors 37  
  correlation matrix 37  
  covariance matrix 37  
  error reduction 37  
  storing intermediate results 37  
Eigenvalues 23  
frequencies after degeneracy 24  
frequencies before degeneracy 24  
linear combinations 23  
of subspectra 23  
of symmetry data 23  
size of submatrices 23  
transitions and energy levels 24

## **P**

PSE window 16

## **R**

resonance frequency 6, 17

## **run**

  generate lineshape 26  
  generate subspectra 26  
  iteration 8, 37  
  simulation 8, 20, 24  
  simulation sequence 25

## **S**

scalar coupling constant 6, 18  
  for aromatic protons of Aspirine 41, 43  
  for epoxy-CH<sub>2</sub>-group of Allylether 53  
  of epoxy-CH"-group of Allylether 52  
simulation of Allyl-glycidyl-ether 48



- 
- experimental spectrum and simulation 55
  - simulation of Asperine spectrum 40
  - smoothing parameters
    - broadening 35
    - iterate baseline 36
    - number of cycles 35
  - spin value 16
  - Spin-System
    - edit 8
    - Edit-Spinsystem window 14
    - export 10
    - load 38
    - open 7, 12
    - store 38
  - start DAISY
    - DAISY command 6
    - from menu bar 6
  - statistical weight 26
  - suppression of transition
    - disable option 17
    - frequency limits 17
    - isolation of overlaid multiplet 17
  - symmetry group 15
    - creation of permutation operators 77
    - Schoenfliess permutation 75
    - Schoenfliess point groups 76

## **T**

- terminate DAISY
  - no save of changes 10, 38
  - save changes 10, 38
- toolbar 7

## **X**

- X-Approximation 16

